



DIGITAL X-RAY CHIROPRACTIC SOLUTIONS

Revolutionary Imaging Systems




2020imaging.net/chiro



Wireless, Fixed, & Hybrid DR



MOMENTUM 14X17 WIRELESS FLATPANEL



- Wireless Interface
- Rapid Charge
- Cesium Scintillator
- Lithium Ion Capacitor
- Featherlight 5.5 lbs
- Cassette Sized
- 16-bit Images

HYBRID 17X17 FLATPANEL



- Wireless or Tethered Interface
- Cesium Scintillator
- Lithium Ion Battery
- Lightweight 9.5 lbs
- Cassette Sized
- 16-bit Images

WIRELESS 14X17 FLATPANEL

- Wireless Interface
- Cesium Scintillator
- Lithium Ion Battery
- Lightweight 8 lbs
- Cassette Sized
- 16-bit Images



TETHERED 17X17 DETECTOR

- Tethered Interface
- Cesium Scintillator
- 16-bit Images
- Cassette Sized



Retrofit and New X-Ray



* RETROFIT ANY X-RAY

Wireless & Fixed DR Compatible



*upper cervical configuration



VERSITILE STRAIGHT-ARM

- Superior Image Quality
- 17x17 Fixed DR compatible
- Space-Saving Design
- Full Range of Motion
- Motorized Vertical Movement
- Detents on 30°, 60°, 90°, & 120°
- H.T. Cables
- 125Kvp, 300k HU, 0.6-1.2mm FS

2020imaging.net/straight-arm

Opal-CHIRO Software

ACQUISITION INTERFACE



- Under 6 second ⌚ Processing
- Preset X-ray Exams
- Anatomical Program (APR)
- Advanced Image Processing
- Touchscreen Enabled
- Background Submission
- Adult & Pediatrics Available
- Fully Customizable
- Image Adjustment Tools
- User Friendly Interface

ADVANCED CHIROPRACTIC TOOLSET

- Advanced Orthogonal Analysis
- Gonstead Measurements
- Blair
- George's Lines
- Ilium Analysis
- Line Drawing
- Measuring
- Angles (Regular, Cobb)
- Cervical Curve
- Lumbar Curve
- Orthospinology
- Full Spine Stitching
- Spine Labeling
- Penning Analysis

Examples of Full Spine Stitching & Gonstead Measurement a



Multi-funtion Studylist

- EHR Integration Ready*
- Burn/Import CDs
- Custom Search/Query
- Customized User Settings



MOBILE VIEWING WITH M-OPAL[†]

iPad®, iPhone®, Android™, Blackberry®, any mobile or smart device!



- Lightweight & Fast
- Patient List & Viewer
- Simple Annotations Enabled
- Tap, Drag, and Pinch
- Patient Filter & Search



CLOUD-BACKUP OPTIONAL

- Live Synchronization
- Fully Monitored
- 448-bit encryption (traffic)
- Fully Customizable
- Enterprise Grade
- Recover Deleted Files



Advanced
Superior
Efficient
Reliable
Affordable

2020imaging.net/chiro

