




DIGITAL X-RAY PODIATRY SOLUTIONS

Revolutionary Imaging Systems



2020imaging.net/podiatry

Wireless and Tethered DR

NEW WIRELESS MOMENTUM DR PANEL

2020imaging.net/momentum



MFG BY
 KONICA MINOLTA

UNIVERSAL DR 1012V PANEL

2020imaging.net/1012v



 Wireless Interface	 Lithium Ion Capacitor (7-yr warranty)
 Featherlight 3.30 lbs	 Rapid Charge
 Cesium Scintillator	 16-bit Images
 100 micron 5 lp/mm	 10x12" Cassette Sized

 Tethered Interface	 Detachable & Durable Cable
 Lightweight 5.5 lbs	 10x12" Cassette Sized
 Cesium Scintillator	 16-bit Images
 125 micron 4 lp/mm	

New and Retrofit X-Ray

NEW PXS710D X-RAY

-  High Frequency Output*
-  Dose Efficient/Ext Tube Life
-  Superb Image Quality
-  Bi-directional Tube Movement
-  Integrated Cable Management
-  Enhanced Patient Support
-  UL Certified 450lbs weight limit
-  Qualifies for ADA Tax Credit
-  UL & CUL Listed
-  Health Canada Approved

New X-Cel X-Ray models also available

RETROFIT ANY X-RAY

Wireless & Tethered DR Compatible
All options qualify for ADA Tax Credit



Example of older X-cel MB-700 modified with new 3.5'h SD Base and 1012V panel.

* Compared to standard frequency (single phase) x-ray.

Hospital-Grade Software

ACQUISITION INTERFACE



- Under 6 second ⌚ Processing
- Background Submission
- Preset X-ray Exams
- Adult & Pediatric Settings
- Annotational Program (APR)
- Fully Customizable
- Advanced Image Processing
- Image Enhancement Tools
- Touchscreen Enabled
- User Friendly Interface

ADVANCED OPAL-RAD VIEWER



Multi-function Studylist ⬇️ Advanced Viewer ➡️

- EHR Integration Ready**
- Podiatric Toolset Available†
- Burn/Import CDs
- Fast, User Friendly Interface
- Custom Search/Query
- Image Fine Tuning
- Customized User Settings
- Exam Room Viewing



MOBILE VIEWING WITH M-OPAL†

iPad®, iPhone®, Android™, Blackberry®, any mobile or smart device!



- Lightweight & Fast
- Patient List & Viewer
- Simple Annotations Enabled
- Tap, Drag, and Pinch
- Patient Filter & Search



CLOUD-BACKUP OPTIONAL

- Live Synchronization
- 448-bit Encryption (traffic)
- Enterprise Grade
- Fully Monitored
- Fully Customizable
- Recover Deleted Files



Advanced
Superior
Efficient
Reliable
Affordable

2020imaging.net/podiatry

