



DIGITAL X-RAY

VETERINARY SOLUTIONS

Revolutionary Imaging Systems



2020imaging.net/vet



Wireless, Fixed, & Hybrid DR



MOMENTUM 14X17 WIRELESS FLATPANEL



- ✓ 3 Year Warranty
- 14"x17" Imaging Area
- Cesium Scintillator
- 100 μ m | 5 lp/mm
- Wireless connectivity
- Cassette Sized
- 7-year lithium-ion capacitor

MOMENTUM-DR

HYBRID 17X17 FLATPANEL



- ✓ 5 Year Warranty
- 17"x17" Imaging Area
- Cesium Scintillator
- 139 μ m | 3.6 lp/mm
- Wireless and tethered connectivity
- Cassette Sized

UNIVERSAL DR

WIRELESS 14X17 FLATPANEL

- ✓ 5 Year Warranty
- 14"x17" Imaging Area
- Cesium Scintillator
- 150 μ m | 3.3 lp/mm
- Wireless connectivity
- Cassette Sized



UNIVERSAL DR

TETHERED 17X17 FLATPANEL

- ✓ 5 Year Warranty
- 17"x17" Imaging Area
- Cesium Scintillator
- 139 μ m | 3.6 lp/mm
- Fixed Tethered Interface
- Cassette Sized



UNIVERSAL DR



Retrofit and New X-Ray



* RETROFIT ANY X-RAY

Wireless & Fixed DR Compatible



⬇ VETSMART™ SYSTEM

- Fixed 40" SID
- Toshiba x-ray tube 140 HU 1-2mm focal spot
- kVP range: 40-150kV
- Single Phase 208-230VAC
- HF 32kW capacitor assisted 400mA generator
- 15" Touch Screen Operator Console
- Six exclusive table base/table top designs
- Push and open tray to access flatpanel
- Exclusive Push and Open Front panel latch mechanism for easy generator access

2020imaging.net/vetsmart

Veterinary Specific Software

ACQUISITION INTERFACE



- | | |
|-------------------------------|---------------------------|
| • Under 6 second ⌚ Processing | • Background Submission |
| • Preset X-ray Exams | • Animal Specific |
| • Anatomical Program (APR) | • Fully Customizable |
| • Advanced Image Processing | • Image Adjustment Tools |
| • Touchscreen Enabled | • User Friendly Interface |

Multi-funtion Studylist ⬇️ Advanced Viewer ➡️

- | | |
|----------------------------|------------------------------------|
| • PM Software Integration* | • Vertibral Heart Score Tool incl. |
| • Burn/Import CDs | • Fast, User Friendly Interface |
| • Custom Search/Query | • Image Fine Tuning |
| • Customized User Settings | • Exam Room Viewing |



MOBILE VIEWING WITH M-OPAL[†]

iPad®, iPhone®, Android™, Blackberry®, any mobile or smart device!



- | |
|------------------------------|
| • Lightweight & Fast |
| • Patient List & Viewer |
| • Simple Annotations Enabled |
| • Tap, Drag, and Pinch |
| • Patient Filter & Search |



CLOUD-BACKUP OPTIONAL

- | | | |
|------------------------|--------------------------------|-------------------------|
| • Live Synchronization | • 448-bit Encryption (traffic) | • Enterprise Grade |
| • Fully Monitored | • Fully Customizable | • Recover Deleted Files |



Advanced
Superior
Efficient
Reliable
Affordable

2020imaging.net/vet

